



Wedge

Wheelchair Cushion



- ☞ Up to stage I wound
- ☞ Helps reduce mild sacral sliding and leaning
- ☞ Moisture-resistant, 4-way stretch cover included



Product Description

HCPC E2601

The Journey™ Wedge cushion is designed to offer comfort and support while reducing sacral sliding. Visco-elastic memory foam conforms to bony prominences while wedge design keeps an individual's pelvis in neutral pelvic alignment. The Wedge reduces extensor tone by increasing hip and knee flexion. Includes a 4-way stretch, incontinent-resistant cover with non-skid bottom. Extra covers available.

Model and Product Specifications

Model #	Size (WIDTH × DEPTH × HEIGHT)	Weight Limit	
SW4-16×16	16" × 16" × 4" front to 2" back	300 lbs.	
SW4-16×18	16" × 18" × 4" front to 2" back		
SW4-18×16	18" × 16" × 4" front to 2" back		
SW4-18×18	18" × 18" × 4" front to 2" back		
SW4-20×16	20" × 16" × 4" front to 2" back		
SW4-20×18	20" × 18" × 4" front to 2" back		

*Other sizes available

Specifications are subject to change without notice.

Copyright © 2010, Keen Mobility Company. All rights reserved. Rev 100122. All trademarks are property of their respective owner.