



Solid Seat Insert

Wheelchair Accessories



- ☞ Creates a level base of support for slung seats
- ☞ Fully covered in non-skid material
- ☞ **MUST** be used with a wheelchair cushion



Product Description

The Solid Seat Insert is designed to offer a flat surface for wheelchairs with slinging or hammocking seats. Fully covered in non-skid material, this insert will stay in place. Use this insert with a Journey™ cushion to help eliminate internal rotation of the hips. Insert is properly labeled to help avoid putting an individual directly onto the insert. Insert must be used with a wheelchair cushion. We recommend ordering one size smaller than seat pan.

Model and Product Specifications

Do NOT PLACE INDIVIDUAL DIRECTLY ON INSERT

Weight Limit: 600lbs.

Model #	Size (WIDTH × DEPTH)
SSI-16×16	16" × 16"
SSI-16×18	16" × 18"
SSI-18×16	18" × 16"
SSI-18×18	18" × 18"
SSI-20×18	20" × 18"