

# AntiThrust


Wheelchair Cushion



- ☞ Up to stage II wound
- ☞ Counteracts forward sliding and mild sacral sitting
- ☞ Moisture-resistant, 4-way stretch cover included




Visco-Elastic Memory Foam  
 High-Resiliency Pressure-Redistributing Foam



(Back View of Cushion)

Ischial Well



## Product Description

HCPC E2601 (3" Only)

The AntiThrust cushion is designed to prevent ischials from thrusting forward. The ischial well is filled with visco-elastic memory foam, which comfortably reduces movement in mild sliders. The unique flat-top design facilitates foot propulsion of the wheelchair while minimizing the risk for hamstring pain or tightness. Includes a 4-way stretch, incontinent-resistant cover with non-skid bottom. Extra covers available.

## Model and Sizes Available

AntiThrust 3" HEIGHT			AntiThrust 4" HEIGHT		
Model #	Size (WIDTH × DEPTH)	Weight Limit	Model #	Size (WIDTH × DEPTH)	Weight Limit
ST3-16×16	16" × 16"	275 lbs.	ST4-16×16	16" × 16"	300 lbs.
ST3-16×18	16" × 18"		ST4-16×18	16" × 18"	
ST3-18×16	18" × 16"		ST4-18×16	18" × 16"	
ST3-18×18	18" × 18"		ST4-18×18	18" × 18"	
ST3-20×18	20" × 18"		ST4-20×18	20" × 18"	

Ischial Well  
Back      Front

Ischial Well  
Back      Front

\*Other sizes available

Specifications are subject to change without notice.

Copyright © 2010, Keen Mobility Company. All rights reserved. Rev 100122. All trademarks are property of their respective owner.